

3D Slicer

Welcome

 Add Data

 Add DICOM Data

 Install Extensions

 Download Sample Data

 Customize Slicer

 Explore Added Data

Feedback

 Share your stories with us on the [Slicer forum](#) and let us know about how 3D Slicer has enabled your research.

We are always interested in improving 3D Slicer. To tell us about your problem or submit a bug report, click [Bug](#).

► About

► Documentation & Tutorials

► Updates

► Acknowledgment

Data Probe

☒ Show Zoomed Slice

L
F
B



3D Slicer

Welcome

 Add Data

 Add DICOM Data

 Install Extensions

 Download Sample Data

 Customize Slicer

 Explore Added Data

Feedback

 Share your stories with the [3D Slicer forum](#) and let us know about how 3D Slicer has enabled your research.

We are always interested in improving 3D Slicer, to tell us about your problem or submit a bug report, open **Help -> Report a Bug**.

▸ About

▸ Documentation & Tutorials

▸ Updates

▸ Acknowledgment

Data Probe

☒ Show Zoomed Slice

L
F
B

R

S: 0.0000mm

G

A: 0.0000mm

1

S

R

P

L

I

L: 0.0000mm

1. Add the Data.

3D Slicer

Welcome

Add Data Add DICOM Data

Install Extensions Download Sample Data

Customize Slicer Explore Added Data

Feedback

Share your stories with us on the [Slicer forum](#) to let us know about how 3D Slicer has enabled your research.

We are always interested in improving 3D Slicer, to tell us your problem or submit a bug report, open **Help -> Report Bug**.

► About

► Documentation & Tutorials

► Updates

► Acknowledgment

Data Probe

☒ Show Zoomed Slice

L
F
B



4 data into the scene

Choose Directory to Add Choose File(s) to Add ☐ Show O


File	Description
✓ /Users/chriswith/Desktop/2120-5-2024-12-18-12-35-39_1.nii	Volume



Reset Cancel **Ok**



2. Add the File with Volume set as the Description.





3. Switch the Screen Setup to one slice. Open the Segment Editor and Add 2 Segments.











6. Go to Grow From Seeds and Initialize. (It may take a while, so be patient)









8. Once the segments are properly formatted, go back to Grow From Seeds and click Apply.




9. Apply Smoothing at 0.1 – 0.2 mm. (Can click Apply multiple times for further smoothing)







11. Go back to the four-screen menu. Click Erase Inside, Rectangle, and Unlimited. Hold shift and move your cursor on the axial view image until the coronal image shows the anterior beginning/formation of the TBI. Stop holding shift and cut out the top half of the axial view image.



12. Hold shift again and move your cursor on the axial image so that the coronal Image shows the posterior ending of the TBI. Stop holding shift and cut out the bottom half of the axial image.



13. Go back to the single screen view. Go to the eraser and change the diameter to around 5% and enable sphere brush. Erase the circular lesion segments into the shape that outlines what a brain would naturally look like if the lesion was not present. (Do this for every slide with the circle)

